

Quadtools Custom Order Form

Each Quadtool is custom-made for each person. Measuring the veteran's arm is a critical part of the process to ensure a proper fit. **Please measure to the nearest 1/8th inch.**

Person Ordering _____

Hand/Arm Left Right

Measurements (rounded to the nearest 1/8")




Forearm Length

Measure the distance from the inside of the elbow to the wrist, with the elbow and wrist bent.



Forearm Circumference

Measure around the thickest part of the forearm, near the elbow.




Wrist & Hand Thickness


Print the following pages, cut along the dotted lines, and use the openings to measure wrist or hand thickness. A clamp may also be used.

The best method for taking wrist thickness and hand thickness measurements is to use any kind of clamp. For the wrist, set the opening of the clamp to just slide over the thinnest part of the wrist, just behind the protruding bone on the outside top of the wrist. For the hand thickness, set the clamp opening to just slide over the middle of the palm. Next, just measure the clamp opening.


Please print these pages and cut along the dotted line, then slide your hand or wrist into the openings to determine the proper size.



Please print these pages and cut along the dotted line, then slide your hand or wrist into the openings to determine the proper size.



Please print these pages and cut along the dotted line, then slide your hand or wrist into the openings to determine the proper size.



Please print these pages and cut along the dotted line, then slide your hand or wrist into the openings to determine the proper size.

